

Portable Ultrasonic Cleaner

NHF/NHPC-B Series



NHF-3010



NHPC-B3010S

Model	NHF-3010	NHPC-B2005S(B)	NHPC-B3010S(B)
Capacity	10 Liter	5.7 Liter	10 Liter
Tank Size (WxDxH)	295x240x160	295x160x160	295x240x160
Overall Size (WxDxH)	490x368x297	325x230x320	325x330x320
Temp. Range	0 ~70 °C		
Timer Range	1 ~ 99 min, 4 hour		
Transducer	BLT Type		
Ultrasonic Frequency	68 KHz, 80KHz,100KHz, 132 KHz Order Made		
Input Power	AC220V 50/60 Hz		
Body Material	ABS	SUS-304(S) or SUS-430 Planting(B)	
Control Method	Microcontroller Control Method		
Switching Method	FET Switching		
Ultrasonic Power (Max)	300W	200W	300W
Heater Power	190W	150W	190W
Weight	11 kg	9 kg	11 kg